


Product Datasheet

IgM Antibody [M-1] (APC) (orb1251548)

Description	IgM Antibody [M-1] (APC)
Species/Host	Mouse
Reactivity	Aves, Gallus
Conjugation	APC
Tested Applications	FC, IHC-Fr
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IgM Antibody [M-1] for use in flow cytometry assays.
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	IgM Antibody [M-1] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



Chicken
peripheral
blood
mononuclear
cel...