


Product Datasheet

CD8B Antibody (PE) (orb1251571)

Description	CD8B Antibody (PE)
Species/Host	Mouse
Reactivity	Feline
Conjugation	PE
Tested Applications	EM, FC, ICC, IHC-Fr, IP, WB
Immunogen	Feline thymocytes
Target	CD8B
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD8 Antibody [fCD8] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P79336
NCBI	P79336
Dilution Range	CD8 Antibody [fCD8] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



Feline
peripheral
blood
lymphocytes
were...