



Product Datasheet

CD4 Antibody (PE) (orb1251585)


Description	CD4 Antibody (PE)
Species/Host	Mouse
Reactivity	Porcine
Conjugation	PE
Tested Applications	FC
Immunogen	Fresh dd miniature swine thymocytes
Target	CD4
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD4 Antibody [74-12-4] for use in flow cytometry assays.
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	CD4 Antibody [74-12-4] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



Porcine
peripheral
blood
lymphocytes
wer...



HeLa cells
transfected
with FLAG
tagged ...



HeLa cells
transfected
with FLAG
tagged ...