


Product Datasheet

H-2Kd/H-2Dd Antibody [34-1-2S] (PE) (orb1251605)

Description	H-2Kd/H-2Dd Antibody [34-1-2S] (PE)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	PE
Tested Applications	FC, IHC-Fr, IP
Immunogen	BDF1 mouse splenocytes
Target	H2-D1; H2-K1
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	H-2Kd/H-2Dd Antibody [34-1-2S] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	H-2Kd/H-2Dd Antibody [34-1-2S] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



BALB/c
mouse
splenocytes
were
stained
wi...