



Product Datasheet

Ly6g Antibody (PE) (orb1251617)


Description	Ly6g Antibody (PE)
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE
Tested Applications	FC, FLISA, ICC, IHC-Fr
Target	Ly6g
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Ly-6G/Ly-6C Antibody [RB6-8C5] for use in flow cytometry assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P35461
NCBI	P35461
Dilution Range	Ly-6G/Ly-6C Antibody [RB6-8C5] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



C57BL/6 mouse bone marrow cells were sta...



Frozen thermally burned mouse tissue sec...



Frozen thermally burned mouse tissue sec...