


Product Datasheet

UBASH3A antibody (orb125164)

Description	Goat Polyclonal to UBASH3A. NP_061834.1; NP_001001895.1, PMID: 15159412 and GenelD: 53347.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	STS2 / TULA (aa470-482)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70kDa band observed in Human Spleen lysates (calculated MW of 69.8kDa according to NP_001001895.1). Recommended concentration: 1-3µg/ml.. IHC: In paraffin embedded Human Breast shows perinuclear staining in the epithelial cells of the tubules. Recommended concentration, 3-5µg/ml.
Protein Sequence	DSGIRISSVFASP
Clonality	Polyclonal
MW	74.1; 69.8
Entrez	53347
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...

