


Product Datasheet

Ptpcr Antibody (APC/Cy7) (orb1251659)

Description	Ptpcr Antibody (APC/Cy7)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	APC/Cy7
Tested Applications	FC, ICC, IHC-Fr
Immunogen	SJL mouse thymocytes and splenocytes
Target	Ptpcr
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD45.1 Antibody [A20] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P06800
NCBI	P06800
Dilution Range	CD45.1 Antibody [A20] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



SJL mouse splenocytes were stained with ...