


Product Datasheet

Cd86 Antibody (PE/Cy5) (orb1251684)

Description	Cd86 Antibody (PE/Cy5)
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE/Cy5
Tested Applications	FC, ICC, IHC-Fr, IHC-P
Immunogen	LPS-activated CBA/Ca mouse spleen B cells
Target	Cd86
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD86 Antibody [GL1] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P42082
NCBI	P42082
Dilution Range	CD86 Antibody [GL1] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



DBA/2 mouse mastocytoma cell line P815 a...