



Product Datasheet

MMP9 antibody (orb12517)


Description	Rabbit polyclonal antibody to MMP9
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A synthetic peptide of human MMP9 (NP_004985.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P14780
Entrez	4318
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse li...



Immunohistochemical staining of human li...



Immunohistochemical staining of mouse he...