


Product Datasheet

APOBEC3D antibody (orb125170)

Description	Goat Polyclonal to APOBEC3D. NP_689639.2, PMID: 16648136 and GenelD: 140564.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	APOBEC3D (aa61-74)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 45kDa band observed in lysates of cell line A549 (calculated MW of 46.6kDa according to NP_689639.2). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Spleen. Recommended concentration: 3.75µg/ml.
Protein Sequence	PKRQSNHRQEVYFR
Clonality	Polyclonal
MW	46.6
Entrez	140564
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 lysate using...