



Product Datasheet

Ptpcr Antibody (APC/Cy7) (orb1251708)

Description	Ptpcr Antibody (APC/Cy7)
Species/Host	Rat
Reactivity	Feline, Human, Mouse
Conjugation	APC/Cy7
Tested Applications	FC, ICC, IHC-Fr, WB
Immunogen	Abelson murine leukemia virus-induced pre-B cell lymphoma
Target	Ptpcr
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD45R Antibody [RA3-6B2] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P06800
NCBI	P06800
Dilution Range	CD45R Antibody [RA3-6B2] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



BALB/c mouse splenocytes were stained wi...



BALB/c mouse splenocytes were stained wi...