

Product Datasheet

Ptpcr Antibody (APC/Cy5.5) (orb1251709)

Description	Ptpcr Antibody (APC/Cy5.5)
Species/Host	Rat
Reactivity	Feline, Human, Mouse
Conjugation	APC/Cy5.5
Tested Applications	FC
Immunogen	Abelson murine leukemia virus-induced pre-B cell lymphoma
Target	Ptpcr
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD45R Antibody [RA3-6B2] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P06800
NCBI	P06800
Dilution Range	CD45R Antibody [RA3-6B2] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.