


Product Datasheet

Ptpcr Antibody (PE/Texas Red) (orb1251715)

Description	Ptpcr Antibody (PE/Texas Red)
Species/Host	Rat
Reactivity	Feline, Human, Mouse
Conjugation	PE/Texas Red
Tested Applications	FC
Immunogen	Abelson murine leukemia virus-induced pre-B cell lymphoma
Target	Ptpcr
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD45R Antibody [RA3-6B2] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P06800
NCBI	P06800
Dilution Range	CD45R Antibody [RA3-6B2] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



BALB/c
mouse
splenocytes
were
stained wi...