



Product Datasheet

TPH2 antibody (orb125175)

Description	Goat Polyclonal to TPH2. NP_775489.2, PMID: 16203956 and GenelD: 121278.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	TPH2 (aa16-29)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in lysates of cell line HEK293 (calculated MW of 56.1kDa according to NP_775489.2). Recommended concentration: 0.5-2µg/ml. Primary incubation was 1 hour.. IHC: Cryosections of Human Hypothalamus. Optimal dilution 1:5000.
Protein Sequence	RGFSLDSAVPEEHQ
Clonality	Polyclonal
MW	56.1
Entrez	121278
Dilution Range	ELISA: 1:128000, WB: 0.5-2 µg/ml, IHC-P: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 lysate u...



Immunohistochemical staining of human hy...