

Product Datasheet

Cd19 Antibody (PE/Cy5.5) (orb1251767)

Description	Cd19 Antibody (PE/Cy5.5)
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE/Cy5.5
Tested Applications	ELISA, FC, ICC, WB
Immunogen	Mouse CD19-expressing K562 human erythroleukemia cells
Target	Cd19
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD19 Antibody [6D5] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P25918
NCBI	P25918
Dilution Range	CD19 Antibody [6D5] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.