


Product Datasheet

NMNAT3 antibody (orb125179)

Description	Goat polyclonal antibody to NMNAT3
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NMNAT3 (aa149-162)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 26kDa band observed in Human Kidney lysates (calculated MW of 24.1kDa according to NP_835471.1). Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	QHNIHLAKEPVQNE
Clonality	Polyclonal
MW	24.1; 18.3
Entrez	349565
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Kidney ly...