


Product Datasheet

Cd4 Antibody (PE/Cy5.5) (orb1251807)

Description	Cd4 Antibody (PE/Cy5.5)
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE/Cy5.5
Tested Applications	FC, ICC, IHC-Fr, IHC-P, WB
Immunogen	Mouse CTL clone V4
Target	Cd4
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD4 Antibody [GK1.5] for use in flow cytometry assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P06332
NCBI	P06332
Dilution Range	CD4 Antibody [GK1.5] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



BALB/c
mouse
splenocytes
were
stained wi...