

Product Datasheet

Vipr1 antibody (orb125181)

Description	Goat polyclonal antibody to Vipr1
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Vipr1 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 52kDa band observed in Mouse and Rat Small Intestine lysates (calculated MW of 52.1kDa according to NP_035833.2). Recommended concentration: 1-3µg/ml.
Protein Sequence	QLFSPIHGYNISRN
Clonality	Polyclonal
MW	52.1
Dilution Range	ELISA: 1:1000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse (A) and R...