


## Product Datasheet

### Tnfrsf1b antibody (orb125187)

|                            |                                                                                                                                                                                                                                                |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal antibody to Tnfrsf1b                                                                                                                                                                                                           |
| <b>Species/Host</b>        | Goat                                                                                                                                                                                                                                           |
| <b>Reactivity</b>          | Mouse, Rat                                                                                                                                                                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                   |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                                                                                                                                                                                 |
| <b>Target</b>              | TNF-R-II / TNFR2 (mouse)                                                                                                                                                                                                                       |
| <b>Preservatives</b>       | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                    |
| <b>Storage</b>             | Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                                                     |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                          |
| <b>Application notes</b>   | WB: Approx 55kDa band observed in Mouse and Rat Spinal Cord lysates (calculated MW of 50.3kDa according to Mouse NP_035740.2). Recommended concentration: 0.2-0.6µg/ml. IHC: Paraffin embedded Human Liver. Recommended concentration: 5µg/ml. |
| <b>Protein Sequence</b>    | KHFCNKTSDTV CAD                                                                                                                                                                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                     |
| <b>MW</b>                  | 50.3                                                                                                                                                                                                                                           |
| <b>Dilution Range</b>      | ELISA: 1:8000, WB: 0.2-0.6 µg/ml, IHC-P: 5 µg/ml                                                                                                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                |



Western blot analysis of Rat Spinal Cord...