

Product Datasheet

UCHL1 Antibody (orb1251871)

Description

UCLH1 Antibody

Species/Host

Mouse

Reactivity

Bovine, Canine, Guinea pig, Human, Mouse, Porcine, Rabbit, Rat, Sheep, Zebrafish

Conjugation

Unconjugated

Tested

IF, WB

Applications
Immunogen

Native UchL1/PGP9.5 protein isolated from human brain was used as the immunogen for this anti-UchL1 antibody.

Target

UCHL1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Western blot: 0.5-1 ug/ml Immunofluorescence (1-2 ug/ml) The optimal dilution of the anti-UchL1 antibody for each application should be determined by the researcher.

Isotype

IgG2a, kappa

Clonality

Monoclonal

Uniprot ID

P09936

Dilution Range

Western blot: 0.5-1 ug/ml Immunofluorescence (1-2 ug/ml) The optimal dilution of the anti-UchL1 antibody for each application should be determined by the researcher.

Expiration Date

12 months from date of receipt.