

Product Datasheet

Adam12 antibody (orb125190)

Description	Goat Polyclonal to Adam12. NP_031426.2, PMID: 21173277 and GeneID: 8038.
Species/Host	Goat
Reactivity	Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Adam12 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments gave an approx 37kDa band in Mouse and Rat Skeletal Muscle lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 98.5
Protein Sequence	KKTTMEKLRCVHP
Clonality	Polyclonal
MW	98.5
Entrez	8038
Dilution Range	ELISA: 1:4000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.