


Product Datasheet

ACPP Antibody (orb1251979)

Description	ACPP Antibody	 IHC testing of human prostate carcinoma ...
Species/Host	Mouse	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IF, IHC	
Immunogen	Recombinant protein was used as the immunogen for this PSAP antibody.	
Target	ACPP	
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide	
Form/Appearance	Liquid	
Concentration	0.2 mg/mL	
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.	
Note	For research use only	
Application notes	Flow Cytometry: 0.5-1ug/10 ⁶ cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml IHC (FFPE): 0.5-1 ug/ml for 30 minutes at RT (1) Titering of the PSAP antibody may be required for optimal performance.1. FFPE testing requires sections to be boiled in pH6 10mM citrate buffer for 10-20 minutes, followed by cooling at RT for 20 minutes, prior to staining.	
Isotype	IgG1	
Clonality	Monoclonal	
Uniprot ID	P15309	
Dilution Range	Flow Cytometry: 0.5-1ug/10 ⁶ cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml IHC (FFPE): 0.5-1 ug/ml for 30 minutes at RT (1) Titering of the PSAP antibody may be required for optimal performance.1. FFPE testing requires sections to be boiled in pH6 10mM citrate buffer for 10-20 minutes, followed by cooling at RT for 20 minutes, prior to staining.	
Expiration Date	12 months from date of receipt.	