

Product Datasheet

TGM3 antibody (orb125200)

Description	Goat Polyclonal to TGM3. NP_003236.3, PMID: 20716179 and GeneID: 7053.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Transglutaminase 3 (aa598-610)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 80kDa band observed in Human Tonsil lysates (calculated MW of 76.6kDa according to NP_003236.3). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	TLEVLNEARVRKP
Clonality	Polyclonal
MW	76.6
Entrez	7053
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Tonsil ly...