


Product Datasheet

DAO antibody (orb125204)

Description	Goat polyclonal antibody to DAO
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	D-amino-acid oxidase (aa286-298)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 37kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 39.5kDa according to NP_001908.3). Recommended concentration: 0.03-0.01µg/ml.
Protein Sequence	RPQIRLEREQLRT
Clonality	Polyclonal
MW	39.5
Entrez	1610
Dilution Range	ELISA: 1:16000, WB: 0.03-0.01 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Cerebellu...