


Product Datasheet

FBXO5 Antibody (orb1252065)


Description

FBXO5 Antibody

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	A human partial recombinant N-terminal protein fragment was used as the immunogen for the EMI-1 antibody.
Target	FBXO5
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Western blot: 1-2 ug/mlImmunofluorescence: 0.5-1 ug/mlImmunohistochemistry (FFPE): 0.5-1 ug/ml for 30 minutes at RT (1)Prediluted format : incubate for 30 min at RT (2)Optimal dilution of the EMI-1 antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0 for 10-20 min followed by cooling at RT for 20 minutes2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	Q9UKT4
Dilution Range	Western blot: 1-2 ug/mlImmunofluorescence: 0.5-1 ug/mlImmunohistochemistry (FFPE): 0.5-1 ug/ml for 30 minutes at RT (1)Prediluted format : incubate for 30 min at RT (2)Optimal dilution of the EMI-1 antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0 for 10-20 min followed by cooling at RT for 20 minutes2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.



Western blot testing of human 1) HeLa an...



IHC testing of FFPE human ovarian carcin...