


Product Datasheet

Trib1 antibody (orb125207)

Description	Goat Polyclonal to Trib1. NP_653132.1, PMID: 20610816 and GeneID: 10221.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Trib1 (mouse) (aa304-317)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 42kDa band observed in lysates of cell line NIH3T3 (calculated MW of 41.3kDa according to NP_653132.1). Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	PEHVSPKARCLIRS
Clonality	Polyclonal
MW	41.3
Entrez	10221
Dilution Range	ELISA: 1:8000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH3T3 lysate u...