


Product Datasheet

DEFB107A antibody (orb125211)

Description	Goat Polyclonal to DEFB107A. NP_001032757.2; NP_001035795.1, PMID: 16033865 and GeneID: 245910; 5036...
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	DEFB107A / DEFB107B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the calculated size of 7.85kDa according to NP_001032757.2.. IHC: Paraffin embedded Human Colon. Recommended concentration: 5µg/ml.
Protein Sequence	HRALISKRMEGHCE
Clonality	Polyclonal
MW	7.85
Entrez	503614, 245910
Dilution Range	ELISA: 1:32000, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Co...