


Product Datasheet

TSHR antibody (orb125215)

Description	Goat Polyclonal to TSHR. NP_000360.2; NP_001018046.1, PMID: 20929407 and GeneID: 7253.
Species/Host	Goat
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TSHR (aa101-115)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 28kDa band observed in Human Thyroid Gland lysates (calculated MW of 28.4kDa according to NP_001018046.1). Recommended concentration: 0.5-2µg/ml.
Protein Sequence	SKVTHIEIRNTRNLT
Clonality	Polyclonal
MW	86.8; 28.4; 30.8
Entrez	7253
Dilution Range	ELISA: 1:16000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Thyroid g...