


Product Datasheet


BCL6 Antibody (orb1252154)

Description	BCL6 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Amino acids 256-389 were used as the immunogen for the Bcl6 antibody.
Target	BCL6
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 1-2 ug/mlWestern blot: 0.5-1 ug/mlOptimal dilution of the Bcl6 antibody should be determined by the researcher.
Isotype	IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P41182
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 1-2 ug/mlWestern blot: 0.5-1 ug/mlOptimal dilution of the Bcl6 antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.


Human Protein Microarray Specificity Validation



Analysis of HuProt(TM) microarray contain...



Western blot testing of 1) partial recom...



SDS-PAGE Analysis of Purified, BSA-Free ...