


Product Datasheet

Tgm6 antibody (orb125216)

Description	Goat Polyclonal to Tgm6. NP_945345.2, PMID: 21106500 and GeneID: 343641.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Tgm6 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 75kDa band observed in Rat Lung lysates (calculated MW of 79.3kDa according to Mouse NP_945345.2 and 80.4kDa according to Rat XP_230601.5). Recommended concentration: 1-3µg/ml.
Protein Sequence	DKYVDSYGRTLED
Clonality	Polyclonal
MW	79.3; 54.0
Entrez	343641
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Lung lysate...