



Product Datasheet

EPCAM Antibody (orb1252165)

Description	EPCAM Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	A human recombinant partial protein corresponding to amino acids 77-202 (extracellular domain) was used as the immunogen for the Ep-CAM antibody.
Target	EPCAM
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	ELISA: 2-4 ug/ml; order the BSA free formulation for coating WB: 0.5-1 ug/ml Titration of the Ep-CAM antibody may be required for optimal performance.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P16422
Dilution Range	ELISA: 2-4 ug/ml; order the BSA free formulation for coating WB: 0.5-1 ug/ml Titration of the Ep-CAM antibody may be required for optimal performance.
Expiration Date	12 months from date of receipt.



IHC testing of
FFPE human
colorectal
car...



SDS-PAGE
Analysis of
Purified, BSA-
Free ...