

Product Datasheet

FOXA1 Antibody (orb1252170)

Description

FOXA1 Antibody

Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

A human recombinant partial protein corresponding to amino acids 372-472 was used as the immunogen for this FOXA1 antibody.

Target

FOXA1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1ug/10⁶ cellsIF: 0.5-1 ug/mlWB: 0.5-1 ug/mlIHC (FFPE): 0.5-1 ug/ml for 30 min at RTThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the FOXA1 antibody to be titrated up or down for optimal performance.

Isotype

IgG2b, kappa

Clonality

Monoclonal

Uniprot ID


[P55317](https://www.uniprot.org/entry/P55317)

Dilution Range


Flow Cytometry: 0.5-1ug/10⁶ cellsIF: 0.5-1 ug/mlWB: 0.5-1 ug/mlIHC (FFPE): 0.5-1 ug/ml for 30 min at RTThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the FOXA1 antibody to be titrated up or down for optimal performance.

Expiration Date


12 months from date of receipt.



IHC testing of
FFPE human
prostate
carcin...



IHC testing of
FFPE human
bladder
carcin...



SDS-PAGE
Analysis of
Purified, BSA-
Free ...