

Product Datasheet

CFTR Antibody (orb1252177)

Description

CFTR Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

A partial recombinant protein corresponding to amino acids 258-385 from the human protein was used as the immunogen for this CFTR antibody.

Target

CFTR

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

IF: 0.5-1 ug/ml IHC (FFPE): 1-2 ug/ml for 30 min at RT The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the CFTR antibody to be titrated up or down for optimal performance.

Isotype

IgG2b

Clonality

Monoclonal

Uniprot ID


P13569

Dilution Range


IF: 0.5-1 ug/ml IHC (FFPE): 1-2 ug/ml for 30 min at RT The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the CFTR antibody to be titrated up or down for optimal performance.

Expiration Date

12 months from date of receipt.



IHC testing of
FFPE human
pancreas
with ...



SDS-PAGE
Analysis of
Purified, BSA-
Free ...