


Product Datasheet

NGFR antibody (orb125218)

Description	Goat polyclonal antibody to NGFR
Species/Host	Goat
Reactivity	Bovine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	p75NTR / CD271 (aa202-215)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 75kDa band observed in Rat Heart lysates (calculated MW of 45.2kDa according to Human NP_002498.1 and 45.4kDa according to Rat NP_036742.2. The observed molecular weight corresponds to earlier findings in literature with different antibodies (Nykja
Protein Sequence	PEGSDSTAPSTQEP
Clonality	Polyclonal
MW	45.2
Entrez	4804
Dilution Range	ELISA: 1:4000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Heart lysat...