

Product Datasheet

DOC2A antibody (orb125219)

Description	Goat Polyclonal to DOC2A. NP_003577.2, PMID: 18354201 and GenelD: 8448.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	DOC2A (aa233-246)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments gave an approx 50kDa band in Human, Mouse and Rat Brain lysates after 0.3 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 44kD. IHC: Paraffin embedded Human Brain (Cortex). Recommended concentration: 5µg/ml.
Protein Sequence	SCYLKELEQAEQGQ
Clonality	Polyclonal
MW	44
Entrez	8448
Dilution Range	ELISA: 1:16000, WB: 0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.