



Product Datasheet

VWF Antibody (orb1252190)

Description	VWF Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Amino acids 1815-1939 from the human protein were used as the immunogen for this vWF antibody.
Target	VWF
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml Immunohistochemistry (FFPE): 1-2 ug/ml for 30 minutes at RT The optimal dilution of the vWF antibody for each application should be determined by the researcher.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P04275
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml Immunohistochemistry (FFPE): 1-2 ug/ml for 30 minutes at RT The optimal dilution of the vWF antibody for each application should be determined by the researcher.
Expiration Date	12 months from date of receipt.



IHC testing of FFPE human spleen with vW...



SDS-PAGE Analysis of Purified, BSA-Free ...