


Product Datasheet

NKX2-2 Antibody (orb1252196)

Description	NKX2-2 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Human recombinant protein was used as the immunogen for this NOX4 antibody.
Target	NKX2-2
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	IF: 1-2 ug/ml ELISA: suitable; order BSA-free formulation for coating WB: 0.5-1 ug/ml Titration of the NOX4 antibody may be required for optimal performance.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	Q9NPH5
Dilution Range	IF: 1-2 ug/ml ELISA: suitable; order BSA-free formulation for coating WB: 0.5-1 ug/ml Titration of the NOX4 antibody may be required for optimal performance.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...