




Product Datasheet

ENO2 Antibody (orb1252204)

Description	ENO2 Antibody	
Species/Host	Mouse	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	FC, IF, IHC, WB	
Immunogen	Amino acids 416-433 of human Neuron Specific Enolase were used as the immunogen for this NSE antibody.	IHC testing of FFPE human cerebellum wit...
Target	ENO2	
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide	
Form/Appearance	Liquid	
Concentration	0.2 mg/mL	
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.	IHC testing of FFPE human pheochromocyto...
Note	For research use only	
Application notes	Western blot: 0.5-1 ug/mlFlow Cytometry: 0.5-1 ug/million cells in 0.1mlIF: 1-2 ug/mlIHC (FFPE): 0.1-0.2 ug/ml for 30 min at RTThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the NSE antibody to be titrated up or down for optimal performance.	
Isotype	IgG2b	
Clonality	Monoclonal	
Uniprot ID	P09104	Western blot testing of human 1) Y79, 2)...
Dilution Range	Western blot: 0.5-1 ug/mlFlow Cytometry: 0.5-1 ug/million cells in 0.1mlIF: 1-2 ug/mlIHC (FFPE): 0.1-0.2 ug/ml for 30 min at RTThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the NSE antibody to be titrated up or down for optimal performance.	
Expiration Date	12 months from date of receipt.	