

Product Datasheet

TDP2 Antibody (orb1252207)

Description

TDP2 Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

Human recombinant full length protein was used as the immunogen for this TDP2 antibody.

Target

TDP2

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1ug/10⁶ cells in 0.1mLIF: 0.5-1 ug/mL IHC (FFPE): 0.5-1 ug/ml for 30 min at RTThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the TDP2 antibody to be titrated up or down for optimal performance.

Isotype

IgG2b, kappa

Clonality

Monoclonal

Uniprot ID


[O95551](#)

Dilution Range


Flow Cytometry: 0.5-1ug/10⁶ cells in 0.1mLIF: 0.5-1 ug/mL IHC (FFPE): 0.5-1 ug/ml for 30 min at RTThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the TDP2 antibody to be titrated up or down for optimal performance.

Expiration Date


12 months from date of receipt.



IHC testing of
FFPE human
duodenal
carci...



IHC testing of
FFPE human
colon
carcinom...



IHC testing of
FFPE human
prostate
carci...