



Product Datasheet

EPO Antibody (orb1252225)

Description	EPO Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	A human partial recombinant protein was used as the immunogen for the EPO antibody.
Target	EPO
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	WB: 0.5-1 ug/mlIF: 1-2 ug/mlImmunohistochemistry (FFPE): 4-8 ug/ml for 30 min at RT (1)Flow Cytometry: 0.5-1 ug/million cells in 0.1mlTiterting of the EPO antibody may be required for optimal performance.1. Required HIER: boil tissue sections in 10mM citrate buffer, pH 6, for 10-20 min.
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P01588
Dilution Range	WB: 0.5-1 ug/mlIF: 1-2 ug/mlImmunohistochemistry (FFPE): 4-8 ug/ml for 30 min at RT (1)Flow Cytometry: 0.5-1 ug/million cells in 0.1mlTiterting of the EPO antibody may be required for optimal performance.1. Required HIER: boil tissue sections in 10mM citrate buffer, pH 6, for 10-20 min.
Expiration Date	12 months from date of receipt.



Western blot testing of with EPO antibod...



SDS-PAGE Analysis of Purified, BSA-Free ...