


Product Datasheet

EPO Antibody (orb1252226)

Description	EPO Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	A human partial recombinant protein was used as the immunogen for the Erythropoietin antibody.
Target	EPO
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	IF: 1-2 ug/ml Immunohistochemistry (FFPE): 4-8 ug/ml for 30 min at RT Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Titration of the Erythropoietin antibody may be required for optimal performance.
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P01588
Dilution Range	IF: 1-2 ug/ml Immunohistochemistry (FFPE): 4-8 ug/ml for 30 min at RT Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Titration of the Erythropoietin antibody may be required for optimal performance.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...