


Product Datasheet

TP63 Antibody (orb1252230)

Description	TP63 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	Recombinant partial human p53 protein was used as the immunogen for this p63 antibody.
Target	TP63
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 1-2 ug/ml Immunohistochemistry (FFPE): 1-2 ug/ml for 30 minutes at RT The optimal dilution of the p63 antibody for each application should be determined by the researcher.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H3D4
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 1-2 ug/ml Immunohistochemistry (FFPE): 1-2 ug/ml for 30 minutes at RT The optimal dilution of the p63 antibody for each application should be determined by the researcher.
Expiration Date	12 months from date of receipt.



IHC testing of FFPE human prostate cancer...