


Product Datasheet

EPCAM Antibody (orb1252236)

Description	EPCAM Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A human partial recombinant protein corresponding to amino acids 77-202 (extracellular domain) was used as the immunogen for the EpCAM antibody.
Target	EPCAM
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	ELISA: order the BSA free formulation for coating Titering of the EpCAM antibody may be required for optimal performance.
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P16422
Dilution Range	ELISA: order the BSA free formulation for coating Titering of the EpCAM antibody may be required for optimal performance.
Expiration Date	12 months from date of receipt.



IHC testing of
FFPE human
ovarian
carcin...