


Product Datasheet

LHCGR Antibody (orb1252240)

Description	LHCGR Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	A partial recombinant protein from the N-terminal region of human Luteinizing hormone / Choriogonadotropin receptor was used as the immunogen for the hCG Receptor antibody.
Target	LHCGR
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1mlIF: 1-2 ug/mlWB: 0.5-1 ug/mlOptimal dilution of the hCG Receptor antibody should be determined by the researcher.
Isotype	IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P22888
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1mlIF: 1-2 ug/mlWB: 0.5-1 ug/mlOptimal dilution of the hCG Receptor antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...