



Product Datasheet

TNFSF15 Antibody (orb1252245)


Description	TNFSF15 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	A full length human recombinant protein was used as the immunogen for this TL1A antibody.
Target	TNFSF15
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1ug/10 ⁶ cellsIF: 1-2 ug/mlWestern blot: 1-2 ug/mlIHC (FFPE): 1-2 ug/mlTitration of the TL1A antibody may be required for optimal performance.
Isotype	IgG1, kappa
Clonality	Monoclonal
Uniprot ID	O95150
Dilution Range	Flow Cytometry: 0.5-1ug/10 ⁶ cellsIF: 1-2 ug/mlWestern blot: 1-2 ug/mlIHC (FFPE): 1-2 ug/mlTitration of the TL1A antibody may be required for optimal performance.
Expiration Date	12 months from date of receipt.



IHC testing of
FFPE human
colon
carcinom...



IHC testing of
FFPE human
parathyroid
ma...



SDS-PAGE
Analysis of
Purified, BSA-
Free ...