


Product Datasheet

GP2 Antibody (orb1252256)

Description	GP2 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	A partial human recombinant protein corresponding to amino acids 35-179 was used as the immunogen for the GP2 antibody.
Target	GP2
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1ug/10e6 cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml IHC (FFPE): 1-2 ug/ml for 30 min at RT Optimal dilution of the GP2 antibody should be determined by the researcher.
Isotype	IgG2b, lambda
Clonality	Monoclonal
Uniprot ID	P55259
Dilution Range	Flow Cytometry: 0.5-1ug/10e6 cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml IHC (FFPE): 1-2 ug/ml for 30 min at RT Optimal dilution of the GP2 antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.



IHC testing of FFPE human pancreas with ...