


Product Datasheet

CDH2 Antibody (orb1252267)

Description	CDH2 Antibody	
Species/Host	Mouse	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	FC, IF, IHC-P, WB	
Immunogen	A human partial recombinant protein corresponding to the intracellular portion of CDH2 was used as the immunogen for this Cadherin 2 antibody.	IHC testing of FFPE mouse heart with Cad...
Target	CDH2	
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide	
Form/Appearance	Liquid	
Concentration	0.2 mg/mL	
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.	
Note	For research use only	
Application notes	Flow Cytometry: 0.5-1ug/10 ⁶ cells in 0.1mlIF: 1-2 ug/mlWB: 0.5-1 ug/mlIHC (FFPE): 0.5-1 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Cadherin 2 antibody to be titrated up or down for optimal performance.	
Isotype	IgG1, kappa	
Clonality	Monoclonal	
Uniprot ID	P19022	
Dilution Range	Flow Cytometry: 0.5-1ug/10 ⁶ cells in 0.1mlIF: 1-2 ug/mlWB: 0.5-1 ug/mlIHC (FFPE): 0.5-1 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Cadherin 2 antibody to be titrated up or down for optimal performance.	
Expiration Date	12 months from date of receipt.	