

## Product Datasheet

### ACPP Antibody (orb1252269)

## Description

ACPP Antibody

|                            |                                                                                                                                                                                                  |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species/Host</b>        | Mouse                                                                                                                                                                                            |
| <b>Reactivity</b>          | Human                                                                                                                                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                     |
| <b>Tested Applications</b> | FC, IF, IHC                                                                                                                                                                                      |
| <b>Immunogen</b>           | PSAP purified from human seminal plasma was used as the immunogen for this PSAP antibody.                                                                                                        |
| <b>Target</b>              | ACPP                                                                                                                                                                                             |
| <b>Preservatives</b>       | PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide                                                                                                                                               |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                           |
| <b>Concentration</b>       | 0.2 mg/mL                                                                                                                                                                                        |
| <b>Storage</b>             | Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.                                                                                                                                             |
| <b>Note</b>                | For research use only                                                                                                                                                                            |
| <b>Application notes</b>   | Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 0.5-1 ug/mlIHC (FFPE): 0.5-1 ug/ml for 30 minutes at RTTitering of the PSAP antibody may be required for optimal performance. |
| <b>Isotype</b>             | IgG1, kappa                                                                                                                                                                                      |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                                       |
| <b>Uniprot ID</b>          | <b>P15309</b>                                                                                                                                                                                    |
| <b>Dilution Range</b>      | Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 0.5-1 ug/mlIHC (FFPE): 0.5-1 ug/ml for 30 minutes at RTTitering of the PSAP antibody may be required for optimal performance. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                  |