

Product Datasheet

GJB1 Antibody (orb1252277)

Description

GJB1 Antibody

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-Fr, WB

Immunogen

Rat junctional complexes were used as the immunogen for the Connexin 32 antibody. The epitope is known to be located in the cytoplasmic loop (Ref 1).

Target

GJB1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 1-2 ug/ml
WB: 0.5-1 ug/ml
IHC (Frozen): 0.5-1 ug/ml for 30 min at RT
Optimal dilution of the Connexin 32 antibody should be determined by the researcher.

Isotype

IgG

Clonality

Monoclonal

Uniprot ID


P08034

Dilution Range

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 1-2 ug/ml
WB: 0.5-1 ug/ml
IHC (Frozen): 0.5-1 ug/ml for 30 min at RT
Optimal dilution of the Connexin 32 antibody should be determined by the researcher.

Expiration Date

12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...