


Product Datasheet

CLTCL1, CLTC Antibody (orb1252284)

Description	CLTCL1, CLTC Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	A recombinant N-terminal protein fragment from human Clathrin heavy chain was used as the immunogen for the Clathrin Heavy Chain antibody.
Target	CLTCL1, CLTC
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 1-2 ug/ml WB: 0.5-1 ug/ml Optimal dilution of the Clathrin Heavy Chain antibody should be determined by the researcher.
Isotype	IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P53675, Q00610
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 1-2 ug/ml WB: 0.5-1 ug/ml Optimal dilution of the Clathrin Heavy Chain antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...